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GENERAL CHARACTERISTICS
Renamel® Microfill is the only composite that truly mimics the enamel surface. Where enamel is to be replicated, Renamel Microfill is the composite resin of choice. Renamel Microfill has been and continues to be the number one rated microfill for more than two decades. When applied, sculpted, finished and polished accurately, no other composite resin can match its lifelike, natural appearance. In essence, it is the artistic dentist's dream material. Renamel Microfill is intended for the replacement of tooth substance and is used solely by the dentist in the dental practice.

DESCRIPTION
Renamel Microfill is a light-cured, microfilled, highly polishable dental filling material, according to EN ISO 4049 (Type 1/Class 2/ Group 1). Renamel Microfill cures with an LED or Halogen light with a light intensity more than 500mW/cm² at wave- length range of 400-500 nm. Renamel Microfill has a compressive strength of >320 MPa, flexural strength of > 80 MPa and 2 mm depth of cure.

CREATING AN ENAMEL-LIKE SURFACE
Renamel Microfill, and as such can be polished as highly as glass. Using the appropriate instruments, a dentist can shape and texture the surface to exactly duplicate the surrounding dentition. This step should not be taken lightly, for this is the part that gives you the ultimate in esthetic dentistry. Renamel Microfill has the best refractive and reflective index of any product in the marketplace. When finished and polished properly, Renamel Microfill will create results that are so beautiful and long-lasting, that both the dentist and the patient will be astounded.

COMPOSITION
Renamel Microfill contains UDMA and 1,4-Butanediol dimethacrylate (28 – 30%). The fillers contain silicon dioxide and prepolymerized composite resin (up to 60% by volume). Initiators, stabilizers, and pigments are additional ingredients (< 1%). The particle size of the fillers range from 0.04 – 0.2 micron.

Renamel Microfill has been in the marketplace for more than twenty- three years, and its handling, color stability, and long term polishability is unsurpassed.

LONG-TERM COLOR STABILITY
Renamel Microfill color stability that is unmatched in the marketplace. Renamel Microfill holds its shade and polish so that you can rest assured that your patient will be happy for many years.

INDICATIONS FOR USE

- Direct anterior veneers
- Top layer of Class I and II restorations
- Class V anterior Class III and IV restorations
- Class V restorations (cervical caries, root erosion, wedge-shaped defects)
- Closing of diastemas
- Repair of composite and ceramic veneers
- Reshape, lengthen and restoratively reposition anterior teeth

CONTRAINDICATIONS
Patient and user: Product should not be used on patients or by users who are known to be allergic to any of the ingredients. Otherwise, the patient group is not restricted.

ADVERSE REACTIONS
Patient and user: Product can cause allergic reactions. User: Observe warnings.

INTERACTIONS
Do not use Renamel Microfill in combination with eugenol containing materials, as these can impede the polymerization process.

NOTICE
During time-consuming procedures, avoid overexposure from the operatory light to protect against premature polymerization of the material. The light curing unit must be maintained and calibrated according to the manufacturer's specifications.

WARNINGS
Contains skin sensitizing substances (acrylates, Skin Sens. 1 > 0.1%). Observe safety data sheet. Applies to the uncured material: May cause an allergic skin reaction. Avoid contact with the uncured material. Wear gloves/protective clothing/eye protection. If on skin: Wash with plenty of water and soap. If skin irritation or rash occurs: Get medical advice/attention. The uncured material is hazardous to the environment. Avoid release to the environment. Always observe the relevant country-specific health and safety regulations and legal requirements when handling hazardous substances. Any deviation from the processing process described here can jeopardize the biocompatibility, functionality and safety of the final product.

DISPENSING COMPOSITE MATERIAL
Syringe delivery: – Extrude composite onto a mixing pad by turning the handle slowly in a clockwise manner. Slice an increment of proper size with a sterile instrument onto a mixing pad. – Immediately replace syringe cap to prevent polymerization. If the extruded material is not used right away, cover the composite with a light protecting box such as our ResinKeeper®. Unit dose delivery – Insert unit dose into dispensing gun and dispense composite material. Dispose of unit dose after use.

PLACEMENT DIRECTIONS
Prior to the application of Renamel Microfill, the use of opaques or microhybrids, or a combination of both may be necessary to achieve excessive color blockout. If replacing dentin, place a micro- hybrid or nanofill composite first for strength and opacity. Polymerize before the addition of Renamel Microfill.

① Prepare tooth, etch and dry according to manufacturer's instructions.

② Place dental adhesive of choice, such as Cosmedent's Complete™, and polymerize. *Other documents to be observed: → [IU Complete](#)™.*

③ Renamel Microfill should be placed in a light-protecting box, such as the ResinKeeper, to prevent prepolymerization by ambient light. Slice an increment of proper size and roll the composite into a ball, a cylinder, or an elliptical shape (convenience form) as needed.

- ④ Place with a clean metal instrument onto the area to be restored and shaped. Cosmedent's BA and BAL (long and short bladed) and G-3 are excellent instruments for placement.
- Note: For instrumentation, titanium-coated composite instruments are desired because they will prevent composite material from sticking to the instrument.
- ⑤ If resin still sticks to the instrument, the instrument may be scratched or dirty. If your room is excessively warm and some stickiness occurs, simply wipe the instrument with a light alcohol swab and this will reduce stickiness.
- ⑥ Take your time and sculpt carefully. Slightly overcontour every restoration, except in the interproximal areas. Sculpt exactly to contour in these areas and it will greatly aid in developing the final contour of your restoration and shorten your finishing and polishing time. Cosmedent's IPC instruments are excellent for interproximal sculpting.
- ⑦ If placing an incisal shade, the use of a titanium coated instrument, such as Cosmedent's Multi-Purpose instrument, will help you create incisal lobes and grooves in your body composite to make room for the application of your incisal shade.
- ⑧ Polymerize each surface for 20–40 seconds (Halogen or LED Lamp) and final cure each restoration on all surfaces for 60 seconds. Do not try to cure to depths over 2 mm. Note: Dentists tend to under-cure composites. You will get best results if you slightly over-cure.
- ⑨ Finishing should be done with a combination of polishing materials such as diamonds, carbide burs and aluminum oxide discs in coarse and medium grits (Mini FlexiDiscs® from Cosmedent).
- ⑩ High-gloss polish can be achieved using aluminum oxide impregnated instruments, such as fine and superfine flexible discs (Mini FlexiDiscs) and FlexiCups®), followed by Enamelize™ polishing paste applied with (Mini FlexiBuffs®), all from Cosmedent. *Other documents to be observed: → [IU Enamelize](#)™.*

For more detailed polishing techniques, please visit Cosmedent's Learning Center at Cosmedent.com/Learning_Center or call a Cosmedent Customer Care Representative at 1.800.621.6729.

COLOR MATCHING
Is there anything more difficult than matching a tooth shade against an existing dentition? The Renamel Restorative System™ is designed to take the guesswork out of shade selection. It recognizes that there are different zones of color and translucency. It understands that color selection should be about achieving the correct enamel shade with predictability. Cosmedent's Renamel Restorative System will help you easily achieve consistently reliable results when selecting a shade.

Color matching is a difficult science, and many people have trouble with color selection.

– Many materials on the market do not match their shade guides. – It is often difficult for an individual to see shade differences. – The color, opacity, translucency and polishability of many materials include some darker geriatric shades for cervical restorations (A5, A6 and C5) and SuperBrite® shades (S81, S82, and S83) to match intensely bleached teeth, or to create the Hollywood smile. This system is also complemented by unique half shades of A1.5, A2.5, and A3.5, which are commonly missing in other systems. The Renamel Restorative System also contains three vividly translucent incisal shades: light, medium and dark.

SHADE SELECTION
When selecting a shade, one should remember that teeth are usually polychromatic, having three zones: a translucent darker incisal, a lighter central body portion and a darker cervical area. Assessing shade on a single tooth using magnification will help you to easily identify these areas. The body shades in The Renamel Restorative System correspond to the central body zone of the tooth surface. Primary shade selection should be taken against the body area of the tooth. A good suggestion is to remove the incisal and gingival portion on each shade tab which may help you achieve a more accurate match to the body area. You can also make your own shade guide out of composite material.

Some additional hints on shade selection are: – Understand how light affects the overall aspect of your shade selection. Assess shade in fluorescent light (where we view one another most of the time) and augment with other types of light, i.e., shade matching light such as Rite-Lite 2™ (from Cosmedent) or viewing in natural light. – Your materials must accurately match an existing shade guide. If not, make your own shade guide. – You can include the patient in shade changing, but never include them in shade matching.

USEFUL TIPS FOR PLACEMENT
① Use of a clean, dry, gloved finger is very effective in rapidly forming the desired surface prior to your sculpting techniques.

- ② The use of titanium coated instruments is desired to prevent composite pull-back.
- ③ Use a Cosmedent BA or BAL (long or short blade) or a G3 for application and general shaping of the material, and the IPTC or IPLC (long or short blade) for finite refinements.
- ④ The use of Cosmedent's Multi-Purpose instrument is very effective for incisal groove and lobe placement, which allows enough room for addition of the incisal material.
- ⑤ Use a Caliper instrument for measurements.
- ⑥ To avoid contamination, wash and dry your gloves thoroughly prior to placement. Do not use powdered gloves.
- ⑦ Make sure your patient removes all makeup prior to treatment.
- ⑧ When applying composite, place convenience form in one piece, then sculpt off excess. This is much easier than placing too little and then adding more.

HYGIENE / DISINFECTION
– Do not reuse a unit dose compule. – Disinfection of the syringe and the unit dose dispenser can be carried out using a composite Multi-Purpose instrument. – All instruments used should be thoroughly washed and autoclaved.

STORAGE INSTRUCTIONS / SHELF LIFE
The composites Renamel Microfill and SuperBrite spontaneously. Store between 39°F – 75°F (4°C – 24°C). Do not use after the expiration date. See imprim/label for expiration date.

IMPORTANT PRODUCT INFORMATION
Federal (U.S.) law restricts this device to sale by or on the order of a dentist.

Note: The product was developed for use in dentistry and must be used according to the instructions for use. Serious incidents must be reported to DeltaMed GmbH and the responsible authority. For the Summary of Safety and Clinical Performance (SSCP) please consult the Eudamed database (<https://ec.europa.eu/tools/eudamed>) or contact sscp@deltamed.de. *Other documents to be observed: → [IU Enamelize](#)™.*

DISPOSAL INSTRUCTIONS
Dispose of contents/ container in accordance with local regulations.

PACKAGING
Renamel Microfill is available in 4 g syringes or 20 Unit Dose compules of 0.2 grams each.

SHADES
A1, A1.5, A2, A2.5, A3, A3.5, A4, A5, A6, B-Zero, B1, B2, B3, B4, C1, C2, C3, C4, C5, D2, D3, D4, Body White, Light Opaque, Light Incisal, Medium Incisal, Dark Incisal, OpaLite.

SuperBrite Shades: SuperBrite 1, SuperBrite 2, SuperBrite 3

Canada Importer: Clinical Research Dental Supplies And Services Inc. 167 Central Ave., 2nd Floor London, Ontario, Canada, N6A 1M6 Establishment License No. 666

Australian Sponsor: Amalgamdent Dental Supplies Ltd 334 Glen Huntly Road Elsterwick, VIC 3185, Australia

Australian Sponsor: Erskine Products Pty Ltd 5–9 Kylie Street, Macksville, NSW 2447, Australia

What are the problems? – Many materials on the market do not match their shade guides. – It is often difficult for an individual to see shade differences. – The color, opacity, translucency and polishability of many materials include some darker geriatric shades for cervical restorations (A5, A6 and C5) and SuperBrite® shades (S81, S82, and S83) to match intensely bleached teeth, or to create the Hollywood smile. This system is also complemented by unique half shades of A1.5, A2.5, and A3.5, which are commonly missing in other systems. The Renamel Restorative System also contains three vividly translucent incisal shades: light, medium and dark.

SHADE SELECTION
When selecting a shade, one should remember that teeth are usually polychromatic, having three zones: a translucent darker incisal, a lighter central body portion and a darker cervical area. Assessing shade on a single tooth using magnification will help you to easily identify these areas. The body shades in The Renamel Restorative System correspond to the central body zone of the tooth surface. Primary shade selection should be taken against the body area of the tooth. A good suggestion is to remove the incisal and gingival portion on each shade tab which may help you achieve a more accurate match to the body area. You can also make your own shade guide out of composite material.

fr
COSMEDENT ® www.cosmedent.com 1.800.621.6729
Renamel ® Microfill COMPOSITE MICROCHARGÉ PHOTOPOLYMÉRISABLE
MODE D'EMPLOI

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