Renamel® Microhybrid | LIGHT-CURED MICROHYBRID COMPOSITE

INSTRUCTIONS FOR USE

DESCRIPTION

Renamel® Microhybrid is a light-cured, microhybrid filled, radiopaque, highly polishable dental filling material, according to EN ISO 4049 (Type 1/Class 2/ Group 1). Renamel Microhybrid cures with an LED or Halogen light at wavelength range of 400 – 500 nm. It is intended for the replacement of tooth substance and is used solely by the dentist in the dental practice. Renamel Microhybrid has a compressive strength of >340 MPa, flexural strength >120 MPa, radiopacity > 1.0 mm and a depth of cure > 1.5 mm.

COMPOSITION
Renamel Microhybrid consists of UDMA, Bis-GMA and 1,4-Butanediol dimethacrylate (25–26%), fillers (containing strontium glass and silicon dioxide (74–75% by weight and 55–56% by volume)), initiators, stabilizers and pigments (< 1%). The particle size of the fillers range from 0,04 – 3 micron with a mean particle size of 0,7 microns.

INDICATIONS FOR USE

– Anterior restorations (Class III and IV)
– Deep Class V restorations (cervical caries, root erosion, wedge-shaped defects)
– Restorations in the posterior region (Class I and II)
– Veneering of discolored teeth
– Closing of diastemas

CONTRAINDICATIONS

Patient and user: Product should not be used on patients or by users who are known to be allergic to any of the ingredients. Otherwise, the patient group is not restricted.

ADVERSE REACTIONS

Patient and user: Product can cause allergic reactions. User: Observe warnings.

INTERACTIONS

Do not use Renamel Microhybrid in combination with eugenol containing materials, as these can impede the polymerization process.

PREPARATION

To remove surface stains, clean the teeth with pumice and water.

SHADE SELECTION

Natural teeth are usually polychromatic, with three zones: a translucent darker incisal area, a lighter central body area and a darker cervical area.

The body shades of Renamel Microhybrid correspond to the central body zone of the tooth surface. Primary shade selection should be taken against the body area of the tooth. Assessing the shade on a clean, moist tooth using magnification will help to easily identify these areas.

Notice: All Renamel Composites match the Vita classical A1–D4® Shade Guide. The use shade guide must exactly match the Renamel Microhybrid shades. If in doubt, consider making an own composite shade guide.

ISOLATION

To ensure a clean working field, preferably use a rubber dam. Alternatively, cotton rolls plus an evacuator can also be used.

DISPENSING COMPOSITE

Syringe delivery: Extrude composite onto a mixing pad by turning the handle slowly in a clockwise manner. Slice an increment of proper size with a sterile instrument. Immediately replace syringe cap to prevent polymerization.

Unit Dose Delivery: Insert unit dose into dispensing gun and dispense composite material. Dispose of unit dose after use.

If the extruded material is not used right away, cover the composite with a light protecting box (e.g. Cosmedent's Resin Keeper).

ANTERIOR RESTORATIONS: CAVITY PREPARATION

① Conduct preparation according to standard preparation rules for all Class III, IV and V restorations. Place mylar strips if needed.

② Cleanse, etch and dry. Place a dental adhesive of choice, such as Cosmedent's Complete™, and polymerize.
Other documents to be observed: → IFU Complete™.

ANTERIOR RESTORATIONS: PLACEMENT AND CURING

Avoid intense light in the working field!

① In case of tetracycline stains, unsightly non-vital teeth or teeth with metal posts and cores, it may be necessary to use of Cosmedent's Creative Color® Pink Opaque prior to the application of Renamel Microhybrid.
Other documents to be observed: → IFU Creative Color® Pink Opaque.

② Place the composite in incremental layers of about 1 mm, not exceeding 2 mm. Slightly over-contour the restoration – however in the interproximal areas sculpt exactly to contour, which will shorten finishing and polishing time. Sculpt and contour with appropriate instruments. Using Cosmedent's titanium coated instrument helps to prevent composite pull-back and sticking. **Polymerize each layer for 20 – 40 seconds** (Halogen or LED light).

③ Before placing the top layer of composite, your restoration can be further characterized using Cosmedent's Creative Color® Opaquers and Creative Color® Tints.
Other documents to be observed: → IFUs Creative Color® Opaquers and Creative Color® Tints.

④ If ultimate esthetics are desired, overlay with Renamel® Microfill™ in the same shade. The polish of Renamel Microfill will last even longer than the polish of Renamel Microhybrid. Sculpt and contour with appropriate instruments. Using Cosmedent's titanium coated instruments helps to prevent composite pull-back and sticking.
Other documents to be observed: → IFU Renamel® Microfill™.

⑤ **Final light-cure each surface for 60 seconds.**

POSTERIOR RESTORATIONS (CLASS I): CAVITY PREPARATION

① Conduct preparation according to standard preparation rules for Class I restorations.

② Cover very deep areas close to the pulp with a calcium hydroxide material and subsequently use a pressure-resistant cement (e.g. glass ionomer cement). For any used materials, always observe the manufacturer's instructions for use. Do not cover other cavity walls, since they can be used to support the bond with an enamel/dentin adhesive.

③ Cleanse, etch and dry. Place a dental adhesive of choice, such as Cosmedent's Complete™, according to the manufacturer's instructions for use and polymerize.
Other documents to be observed: → IFU Complete™.

POSTERIOR RESTORATIONS (CLASS II): PLACEMENT AND CURING

Avoid intense light in the working field!
① Place Renamel Microhybrid, adding enough composite to cover the entire base of the preparation to a depth about 1 mm and not exceeding 2 mm. Sculpt and contour with appropriate instruments. Using Cosmedent's titanium coated instruments helps to prevent composite pull-back and sticking. **Polymerize each surface for 20 – 40 seconds** (Halogen or LED light).

② Continue to build-up in incremental layers and slightly over-contour the occlusal surfaces.

③ Before placing the top layer of composite, your restoration can be further characterized using Cosmedent's Creative Color® Opaquers and Creative Color® Tints.
Other documents to be observed: → IFUs Creative Color® Opaquers and Creative Color® Tints.

④ **Final light-cure each surface for 60 seconds.**

POSTERIOR RESTORATIONS (CLASS II): CAVITY PREPARATION

① Conduct preparation according to standard preparation rules for Class II restorations.

② Cover very deep areas close to the pulp with a calcium hydroxide material and subsequently use a pressure-resistant cement (e.g. glass ionomer cement). For any used materials, always observe the manufacturer's instructions for use. Do not cover other cavity walls, since they can be used to support the bond with an enamel/dentin adhesive.

③ Cleanse, etch and dry. Place a dental adhesive of choice, such as Cosmedent's Complete™, according to the manufacturer's instructions for use and polymerize.
Other documents to be observed: → IFU Complete™.

④ Place a matrix band, such as Cosmedent's Full Contact or Ultrathin Matrix band, wedge thoroughly and make sure that the gingival margin is completely sealed. Often, it is necessary to wedge from both the buccal and lingual to obtain complete closure.

⑤ Burnish the band using a round ball burnisher and the convex side of a spoon excavator in order to create a contact area and not just a contact point. Use heavy pressure and observe the burnished area closely. You should easily be able to see the contact area you have created.

⑥ Apply a flowable cavity liner such as Cosmedent's hyper-radiopaque Renamel® DE-MARK™ at the gingival seat to differentiate your restoration from tooth structure. Apply a small amount of material and gently tap to place with a small round ball burnisher and cover the gingival seat only. **Light-cure for 40–60 seconds.**
Other documents to be observed: → IFU Renamel® DE-MARK™.

POSTERIOR RESTORATIONS (CLASS II): PLACEMENT AND CURING

Avoid intense light in the working field!

① Place the composite in incremental layers of about 1 mm, not exceeding 2 mm. Make sure that all aspects of the gingival seat are covered and all margins all completely sealed. You can use an IPC Carver to contour the proximal areas to reduce flashing. Once the first layer has been placed, **polymerize for 20–40 seconds** (Halogen or LED light).

② Continue to build-up in incremental layers and slightly over-contour the occlusal surfaces. Note: To prevent the cuspal flexion, it might be useful to develop one cusp at a time.

③ Before placing the top layer of composite, your restoration be further characterized using Cosmedent's Creative Color® Opaquers and Creative Color® Tints.
Other documents to be observed: → IFUs Creative Color® Opaquers and Creative Color® Tints.

④ **Final light-cure each surface for 60 seconds.**

OCCUSAL ADJUSTMENTS

Check occlusion with articulating paper. Pay attention to centric and lateral excursion contacts also. Adjust the occlusal surfaces with fine polishing diamonds or stones.

FINISHING AND POLISHING

Finishing should be done with a combination of polishing materials such as diamonds, carbide burs, and aluminum oxide discs in coarse and medium grits. High-gloss polish can be achieved using fine and superfine flexible discs and cups, followed by Enamelize® polishing paste applied with a suitable buff, such as Cosmedent's Mini FlexiBuff or Felt Flexipoint.

Other documents to be observed: → IFU Enamelize®.

HYGIENE

Do not reuse a unit dose compule. Disinfection of the syringe and the unit dose dispenser can be carried out using a commercially available disinfecting solution. All instruments used should be thoroughly washed and autoclaved.

WARNINGS

Contains skin sensitizing substances (acrylates Skin Sens. 1 > 0,1%), contains reproductive toxic substances (diboron trioxide, Repr. 1B > 0,1%).

Obtain special instructions before use. Observe safety data sheet. Applies to the uncured material: May cause an allergic skin reaction. May damage fertility. May damage the unborn child. Avoid contact with the uncured material. Wear protective gloves/ protective clothing/eye protection.

IF ON SKIN: Wash with plenty of water and soap. IF exposed or concerned: Get medical advice/attention. The uncured material is hazardous to the environment. Avoid release to the environment. Always observe the relevant country-specific health and safety regulations and legal requirements when handling hazardous substances.

Any deviation from the process described here can jeopardize the biocompatibility, functionality and safety of the final product. The light curing unit must be maintained and calibrated according to the manufacturer's specifications.

STORAGE INSTRUCTIONS / SHELF LIFE

Protect from exposure to light. May polymerize spontaneously. Store between 39°F – 75°F (4°C – 24°C). Do not use after the expiration date. See imprint/label for expiration date.

DISPOSAL INSTRUCTIONS

Dispose of contents/container in accordance with local regulations.

DELIVERY UNITS

Renamel Microhybrid is available in 5 g syringes, or packages of 10 Unit Dose compules of 0.25 grams each.

SHADES

A1, A2, A3, A3.5, A4, B-Zero, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4, Body White, Clear Occlusal, Milky White Occlusal, Light Incisal, Medium Incisal, SuperBrite 1, SuperBrite 2, SuperBrite 3.

Rx IMPORTANT PRODUCT INFORMATION

Federal (U.S.) law restricts this device to sale by or on the order of a dentist.

The product was developed for use in dentistry and must be used according to the instructions for use. Serious incidents must be reported to DeltaMed GmbH and the responsible authority. For the Summary of Safety and Clinical Performance (SSCP) please consult

the Eudamed database (<https://ec.europa.eu/tools/eudamed>) or contact sscp@deltamede.de.

Canada Importer:
Clinical Research Dental Supplies And Services Inc.
167 Central Ave., 2nd Floor
London, Ontario, Canada, N6A 1M6
Establishment Licence No. 666

Australian Sponsor:
Amalgadent Dental Pty Ltd
334 Glen Huntly Road
Elsterwick, VIC 3185, Australia

Australian Sponsor:
Erskine Products Pty Ltd
5–9 Kylie Street,
Mackville, NSW 2447, Australia

 Batch number

 Item number

 Medical device

 Shelf life

 Observe instruction for use

 Distributor

 Temperature limit

 Manufacturer

 Keep away from light

 CE MARK

 Contains skin sensitizing substances

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