Insure® / Insure® Lite / Insure® Color Modifier | UNIVERSAL CEMENTATION SYSTEM

INSTRUCTIONS FOR USE

DESCRIPTION

Insure® is a dual cured, microhybrid filled, radiopaque, dental resin cement according to EN ISO 4049 (Type2/Class2/Group1 and Type2/Class3). **Insure resin cement** is available in two viscosities. **Insure resin cement** in the standard kit is formulated in a medium, thixotropic viscosity. **Insure Lite** resin cement is produced in a lower, more flowable viscosity. It is recommended that Insure be used for ceramic laminate veneers and **Insure Lite** be used for ceramic or composite crowns, inlays, onlays, and bridges. Insure cures with light at wavelength range of 400 – 500 nm.

Insure Simulcure Catalyst can be mixed with the colored resin cement in equal parts to produce a complete chemical cure in areas that are less accessible to light-curing. If using the dual-cure catalyst with Insure / Insure Lite the working time is 3 – 4 minutes and setting time is 5 – 6 minutes. **Insure Lite Clear** is available in an **Automix** version.

The shade of **Insure** can be modified with **Insure Color Modifier**.

SHADES

Insure and **Insure Lite** are available in the following shades:

Clear, Yellow Red Light, Yellow Red Universal, Grey Light, White Opaque, Yellow Opaque, and Pink Opaque.

Insure Hi-Chroma Color Modifier is available in White.

RELIABLE & PREDICTABLE

Insure resin cements have over 15 years of reliable performance and top ratings by the major dental materials advisories. This versatile resin cement system allows you to cement in veneers, inlays, onlays, crowns, metal and posts and cores with complete predictability. Insure features easy setting and cleanup, a choice of two viscosities, and a wide color range for predictable chairside adjustments. Insure can be made dual-cure for inlays, onlays and crowns with the addition of **Insure Simulcure Catalyst / Insure Lite Simulcure Catalyst**.

COMPOSITION

Insure contains Bis-GMA, UDMA and 1,4-Butanediol dimethacrylate (30%). The fillers contain barium glass and silicon dioxide (69% by weight and 46% by volume). Initiators and stabilizers are additional ingredients (< 1%). The particle size of the fillers range from 0.04 – 3 micron.

Insure Lite and **Insure Color Modifier** contain Bis-GMA, UDMA and 1,4-Butanediol dimethacrylate (35%). The fillers contain barium glass and silicon dioxide (64% by weight and 41% by volume). Initiators and stabilizers are additional ingredients (< 1%). The particle size of the fillers range from 0.04 – 3 micron.

Insure Simulcure Catalyst / Insure Lite Simulcure Catalyst contain Bis-GMA, UDMA and 1,4-Butanediol dimethacrylate (30% / 35%). The fillers contain barium glass and silicon dioxide (69% / 64% by weight and 46% / 42% by volume). Both contain benzoyl peroxide. The particle size of the fillers range from 0.04 – 3 micron.

INDICATIONS

(Areas of Application)

- Cementation of ceramic veneers
- Cementation of ceramic or composite crowns, inlays, onlays, and bridges

DISPENSING INSURE

Place the disposable tip over the syringe opening and screw it on firmly before dispensing the **Insure** resin cement. This will provide controlled placement into the veneer shell without excess resin cement extruding from the syringe causing material waste. Throw away tip after use and immediately replace the cap on the Insure syringe to avoid pre-polymerization.

EVALUATING FOR MARGINAL FIT

For cementation of ceramic laminate veneers, it is advisable to evaluate each cement laminate veneer for proper proximal fit and interproximal adaptation to adjacent ceramic laminate veneers or natural tooth structure. Any aberrations should be identified, and appropriate measures set up for their elimination.

SELECTING SHADE

We recommend that you use **Prevue® Try-In, Cosmedent's** water soluble matching try-in gels, to select the correct shade of your Insure resin cement. **Prevue Try-In Gels** are especially effective for feldspathic porcelain veneers or any minimal or no-prep veneers. These non-light activated pastes accurately match the polymerized color of our seven resin cements and wash away quickly with water, making shade selection predictable and fast.

PLACEMENT DIRECTIONS

Insure is recommended for ceramic veneers and **Insure Lite** (thin viscosity) is recommended for ceramic or composite crowns, inlays, onlays, bridges and veneers because of the optimum film thickness.

Once the ceramic laminate veneers have been tried in and approved for marginal fit, proximal adaptation, and resin cement color, the following cementation procedure is undertaken:

- ① Isolate the teeth to be veneered with a rubber dam or cotton rolls and an effective evacuation system that will sustain a dry uncontaminated field throughout the duration of the cementation procedure.
- ② Place retraction cord, if needed, along the cervical margins of all teeth to be veneered to permit clear visualization of the cervical finish lines and counteract the potential for soft tissue contamination of any kind.
- ③ Floss each contact to remove any residual provisional veneer resin or resin cement. Relax but do not open the contacts with controlled passage of narrow, Superfine, metal backed **FlexiStrips™** (**Cosmedent**). Remove any residual provisional veneer cement on the facial surface with a fine diamond bur.
- ④ Clean all tooth surfaces to be veneered by initially pumicing all accessible tooth surfaces with pumice. Rinse away all pumice residue and microetch any existing resin restorations to which the veneers will be cemented. Rinse all tooth surfaces to be cemented thoroughly with water and air dry.
- ⑤ If the laboratory has not acid etched and silanated the restorations or if this surface is altered by the dentist, then these steps must be performed prior to resin bonding. If the laboratory has etched the porcelain, skip to step 5c.
 - a) Etch the intaglio of the veneer for 20 seconds with the 5% Porcelain etchant (e.max, Ivoclar Vivadent) or 3 minutes with 9.5% Hydrofluoric acid if Feldspathic porcelain.
 - b) Rinse thoroughly with water and air dry with oil free air to visualize a frosted, etched appearance.
 - c) Apply a silanating agent such as **Insure Silanator (Cosmedent)** for 30 seconds to all etched inner surfaces of the veneers and gently air dry. Place the ceramic laminate veneers in sequential order for cementation in pairs, starting with the central incisors, or in quartets, starting with the incisors.
- ⑥ Place a thin metal strip interproximally on the distal proximal surface of the last prepared tooth to be cemented in this grouping (either a pair of veneers or 4 consecutive veneers starting at the midline), to prevent etching an adjacent tooth not in this cementation grouping. Etch each tooth surface to be veneered with 35%–37% phosphoric acid for approximately 20 seconds. Rinse off all acid with a gentle water spray for 5 seconds and carefully dry the teeth with a clean air spray to a damp, but not desiccated state.
- ⑦ When cementing inlays, onlays, crowns or bridges, mix equal parts of the selected **Insure** base shade in conjunction with the **Insure Simulcure Catalyst** (dual-cure catalyst) to ensure proper polymerization. If you are working with veneers less than 1.5 mm thick, it is preferable to always use light-curable Insure base shade without **Simulcure Catalyst**, otherwise, the material may partially or totally set up before subsequent veneers are ready to be placed. Working time for Insure / Insure Lite is 3 – 4 minutes and setting

time is 5 – 6 minutes.

- ⑧ Load the veneers to be placed with the selected Insure resin cement. Be sure to place the disposable tip over the syringe opening and screw it on firmly before dispensing the Insure resin cement. This will provide controlled placement into the veneer shell without excess resin cement extruding from the syringe causing material waste. Be sure to load each veneer fully to the margin because underloading can cause voids.
- ⑨ Seat the veneer(s) from the incisofacial or facial, using the incisal porcelain as the initial seating guide. Finalize the seating with gentle finger pressure, as the excess Insure resin cement gradually extrudes along the margins. Remove resin cement excess with an **IPC Composite Placement Instrument (Cosmedent)** at the cervical and incisal. Remove all but a narrow string of excess resin cement (0.2 mm wide) at the margin to offset suckback porosity or "crevicing" due to polymerization shrinkage. Spot-cure the veneer with a 2 or 3 mm curing tip on the facial surface away from the margins for 10 seconds to secure it in place. Further refine resin cement removal along the margins as needed with dental floss before light-curing the rest of the veneer. Place glycerin such as **Oxygene™** (**Cosmedent**) along the margins. Light-cure (preferable with 2 lights at once) the veneer for at least 60 seconds facially and 60 seconds incisolingually with a wide diameter curing light. This should be done even if an **Insure / Simulcure** mix was used.
- ⑩ Remove ceramic resin cement flash along the margins with metal-backed **FlexiStrips™** (**Cosmedent**). Bard Parker 12B scalpel blades and **Mini FlexiDiscs™** (**Cosmedent**). If ceramic "ledging" requires adjustment, use microfine diamonds with a copious water spray, **All Ceramic Polishers™** (**Cosmedent**) and polishing paste to renew the glaze.
- ⑪ Check the occlusion in centric and eccentric movements with fine articulating paper. Adjust the occlusion where necessary. Repolish the adjusted surface as suggested in #10.

HYGIENE/DISINFECTION

- Disinfection of the syringe can be carried out using a commercially available disinfecting solution.
- All instruments used should be thoroughly washed and autoclaved.

STORAGE

- Once opened, do not refrigerate
- Keep **Insure** at room temperature (approximately 39 – 75°F; 4 – 24°C)
- Close syringe immediately after use. Exposure to light causes premature polymerization
- Do not use after expiration date
- Apply a silanating agent such as **Insure Silanator (Cosmedent)** for 30 seconds to all etched inner surfaces of the veneers and gently air dry. Place the ceramic laminate veneers in sequential order for cementation in pairs, starting with the central incisors, or in quartets, starting with the incisors.
- Place a thin metal strip interproximally on the distal proximal surface of the last prepared tooth to be cemented in this grouping (either a pair of veneers or 4 consecutive veneers starting at the midline), to prevent etching an adjacent tooth not in this cementation grouping. Etch each tooth surface to be veneered with 35%–37% phosphoric acid for approximately 20 seconds. Rinse off all acid with a gentle water spray for 5 seconds and carefully dry the teeth with a clean air spray to a damp, but not desiccated state.

IMPORTANT PRODUCT INFORMATION

Federal (U.S.) law restricts this device to sale by or on the order of a dentist.

Note: This product was developed for use in the dental sector and should be used in accordance with the instructions for use. The manufacturer does not accept liability for damage caused by its use for any other purpose. The user is also obliged to accept responsibility for checking before use that the material is suitable and can be used for the intended purpose, particularly if this is not listed in the instructions for use.

DISPOSAL INSTRUCTIONS

Dispose of contents/ container in accordance with local regulation.

CONTRAINDICATIONS

Insure is contraindicated if a patient is known to be allergic to any of the ingredients of Insure.

SIDE EFFECTS

May cause irritation in sensitive patients.

INTERACTIONS

Do not use in combination with eugenol containing materials as these may impede the polymerization process.

NOTICE

During time-consuming procedures, avoid overexposure from the operating light to protect against premature polymerization of the material. The light curing unit must be maintained and calibrated according to the manufacturer's specifications. Serious incidents must be reported to DeltaMed GmbH, the responsible importer and the responsible authority.

WARNING

Insure contains 1,4-Butanediol dimethacrylate, Urethane dimethacrylate and Di-benzoyl peroxide (only in **Simulcure Catalyst**). May cause an allergic skin reaction. Causes skin irritation. Causes serious eye irritation. Wear protective gloves. In case of contact with the eyes, rinse immediately with plenty of water and seek medical attention.

PACKAGING

Insure 3 g syringes (includes 3 Nozzle Tips White)
Insure Light 3g syringes (includes 3-17 gauge tips)
Insure Color Modifier (includes 3-19 gauge tips)

Canada Site Licence Holder:
CLINICAL RESEARCH DENTAL
167 Central Ave, Suite 200
London, Ontario, Canada N6A 1M6
Establishment License No. 666

Australia Sponsor:
ERKSINE ORAL CARE
5-9 Kylie St,
Macksville NSW 2447
Australia

Australia Sponsor:
Amalgadent Dental Supplies Pty Ltd
334 Glen Huntly Road
Elsternwick, Victoria 3185 Australia

	Batch number
	Item number
	Medical device
	Observe instruction for use
	Distributor
	Recommended storage temperature
	Manufacturer
	Keep away from light
	Warning

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