Renamel® Flowable Microhybrid

www.cosmedent.com | 800.621.6729

LIGHT-CURED MICROHYBRID COMPOSITE

INSTRUCTIONS FOR USE

GENERAL CHARACTERISTICS
Renamel® Flowable Microhybrid offers a wide shade selection, long-term color stability, and excellent handling properties with a flowable consistency that is great for quick and simple composite repairs and small occlusal repairs of all types. Renamel Flowable Microhybrid is radiopaque to distinguish from tooth structure.

DESCRIPTION
Renamel Flowable Microhybrid is a light cured, microhybrid filled, radiopaque, highly polishable dental filling material according to EN ISO 4049 (Type 1/Class 2/ Group 1). Renamel Flowable Microhybrid cures with an LED or Halogen light with a light intensity more than 500 mW/cm² at a wavelength range of 400 - 500 nm. Renamel Flowable Microhybrid has a compressive strength of > 340 MPa, flexural strength > 100 MPa and a depth of cure > 2 mm.

COMPOSITION
Renamel Flowable Microhybrid contains UDMA, Bis-GMA and 1,4-Butanediol dimethacrylate (32 - 34%). The fillers contain strontium glass and silicon dioxide (65 - 67% by weight and 45 - 46% by volume). Initiators, stabilizers and pigments are additional ingredients (< 1%). The particle size of the fillers range from 0.04 - 3 micron with a mean particle size of 0.7 microns.

INDICATIONS FOR USE
- Class V restorations
- Sealing margins
- Repairing pitted areas in hybrid restorations
- Small occlusal and lingual repairs
- Sealing gingival margins in Class II and Class III restorations
- In deep Class V restorations to block out shine through prior to overlay of Renamel Microfill.

DISPENSING FLOWABLE COMPOSITE
Place the disposable tip over the syringe opening and place it on firmly before dispensing the flowable composite. Throw away tip after use and immediately replace the cap on the syringe to avoid prepolymerization.

PLACEMENT DIRECTIONS
Before treatment, tooth color should be assessed using a standard color guide.
① Keep the teeth moist when assessing the shade. For predictable results, use the same shade as your body composite. Prepare the cavity according to the general standards in composite filling therapy.

② **Apply etching gel to entire cavity, beginning with enamel areas.** Leave the gel for 20 seconds on the tooth structure (total-etch technique). Remove the etching gel and spray the tooth with water for 15 seconds. Make sure no residue remains on the tooth.

③ **Dry cavity with oil-free air.** Re moisten dentin and keep moist. Apply dentin adhesive, e. g. (Complete™ from Cosmedent) according to the Instructions for Use and cure.

④ For deep cavities with only a minimal layer of dentin over the pulp or those which have already had a small area exposed, the pulp should be indirectly or directly capped locally using calcium hydroxide compound. The calcium hydroxide compound should be covered with a glass ionomer cement before the Renamel Flowable Microhybrid is inserted into the cavity.

⑤ **Apply Renamel Flowable Microhybrid in thin layers (max. 2 mm).** Polymerize each layer for 20-40 seconds and final cure each restoration on all surfaces for 60 seconds.

Seal Margins
After placing restoration and sculpting toward margins, polymerize for 40 seconds. Next, gently apply Renamel Flowable Microhybrid at the gingival margin by expressing a slight amount from the syringe along the entire margin and polymerize again, then polish. If a pit is seen at the gingival margin, roughen the area with a diamond bur, etch, wash and dry. Add bonding adhesive, then add a small amount of flowable to the defective area. Polymerize, finish, and polish.

Repair Pitted Areas in Existing Restorations
Roughen the pitted area with a broad round diamond bur, bevel edges, wash, dry, and add bonding adhesive. Syringe flowable, polymerize, finish, and polish.

Occlusal Restorations
Take the smallest bur possible, such as a 330 carbide, remove decay (or use air abrasion or hard tissue laser), wash and dry the area, then etch. Add a primer, place bonding agent, polymerize and syringe flowable into cavity. Make sure the flowable extends slightly past the margin so when you polish and come back to the margin, it will not be under contoured. Remove excess material, polymerize, finish, and polish.

Seal Gingival Seat
To completely seal the gingival seat area prior to placement of your composite, place a thin layer of Renamel Flowable Microhybrid or Renamel De-Mark™ at that gingival seat area and polymerize. This reduces polymerization shrinkage and post-operative sensitivity.

To Start Deep Class V Restorations Where You Need to Block Shine-through
Syringe Renamel Flowable Microhybrid into deep areas of preparation after proper bonding technique and polymerize for 40 seconds. Then, overlay with a final coat of Renamel Microfill.

Temporary Restorations
Slightly reduce and abrade your completed temporary veneer or crown, then wash, clean and dry, place bonding adhesive and polymerize. Add flowable of choice to the surface, polymerize, finish and polish.

Sealing Open Margins - Newly Placed Composite
On a freshly placed composite restoration where the oxygen inhibited layer has not been removed, you can flow material right into the margin and it will adhere to both the dental margin and the restorative material, sealing margin completely. Polymerize, finish and polish as usual.

- Older Polished Composite
On a new or old polished restoration, the margin of the restoration to be sealed must be roughened with a medium diamond bur followed by the total etch technique. Bonding adhesive should be placed and polymerized, then flowable can be added to the margin. Polymerize and polish.

Finishing
Finishing should be done with a combination of polishing materials such as diamonds, carbide burs and aluminum oxide discs in coarse and medium grits (Mini FlexiDiscs® from Cosmedent).

Polishing
High-gloss polish can be achieved using aluminum impregnated instruments fine and superfine flexible discs (Mini FlexiDiscs and FlexiCups®), followed by Enamelize® polishing paste applied with Mini FlexiBuffs®, all from Cosmedent.

For more detailed polishing techniques, please visit Cosmedent's Learning Center at Cosmedent.com/Learning_Center or call a Cosmedent Customer Care Representative at 800-621-6729.

HYGIENE/DISINFECTION
Disinfection of the syringe can be carried out using a commercially available disinfecting solution. All instruments used should be thoroughly washed and autoclaved.

STORAGE
- Once opened, do not refrigerate.
- Keep Renamel Flowable Microhybrid at room temperature (approximately 39 - 75 °F; 4 - 24 °C).
Do not use after expiration date.
- Keep away from sunlight.
- Keep material out of reach of children.

IMPORTANT PRODUCT INFORMATION
Federal (U.S.) law restricts this device to sale by or on the order of a dentist.

Note: This product was developed for use in the dental sector and should be used in accordance with the instructions for use. The manufacturer does not accept liability for damage caused by its use for any other purpose. The user is also obliged to accept responsibility for checking before use that the material is suitable and can be used for the intended purpose, particularly if this is not listed in the instructions for use.

DISPOSAL INSTRUCTIONS
Dispose of contents/ container in accordance with local regulation.

CONTRAINDICATIONS
Renamel Flowable Microhybrid is contraindicated if patient is known to be allergic to any of the ingredients of Renamel Flowable Microhybrid. Renamel Flowable Microhybrid is contraindicated in posterior restorations in contact areas.

SIDE EFFECTS
Can cause irritation in sensitive patients.

INTERACTIONS
Do not use Renamel Flowable Microhybrid in combination with eugenol containing materials, as these can impede the polymerization process.

NOTICE
During time-consuming procedures, avoid overexposure from the operatory light to protect against premature polymerization of the material. The light curing unit must be maintained and calibrated according to the manufacturer's specifications. Serious incidents must be reported to DeltaMed GmbH, the responsible importer and the responsible authority.

WARNING
Renamel Flowable Microhybrid contains 1,4-Butanediol dimethacrylate and Urethane dimethacrylate. May cause an allergic skin reaction. Wear protective gloves. In case of contact with the eyes, rinse immediately with plenty of water and seek medical attention.

PACKAGING
Renamel Flowable Microhybrid is available in 3 g syringes, (includes 3-22-gauge tips).

SHADES
A1, A2, A3, A3.5, A4
B1, B2, B3
C2, C3, C4, D3
Body White, Medium Incisal
Clear Occlusal, Milky White Occlusal

Canada Site Licence Holder:
CLINICAL RESEARCH DENTAL
167 Central Ave, Suite 200
London, Ontario, Canada N6A 1M6
Establishment License No. 666

Australia Sponsor:
Amalgadent Dental Supplies Pty Ltd
334 Glen Hunthly Road
Elsternwick, Victoria 3185 Australia

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Cosmedent® Inc.