








COMPOSITE POLISHING GUIDE



Cosmedent gives you everything you need to polish composites. With our conveniently color-coded easy to use polishers, you can create a smooth and esthetic polish that mimics enamel long-term.

Cosmedent Polishing Tool	Speed Range	Use Notes	Sizes
Coarse FlexiDisc® 	3,000 to 15,000 RPM for 3/8" disc 5,000 to 15,000 RPM for all other sizes	<ul style="list-style-type: none"> • Use light pressure • Gross reduction • Initial contouring • Buffing out imperfections 	<ul style="list-style-type: none"> • Available Regular center sizes: 1/2", 5/8", & 3/4" • Available Mini center sizes: 3/8", 1/2", & 5/8"
Medium FlexiDisc® 	3,000 to 10,000 RPM around gingival margin with 3/8" disc 10,000 to 15,000 RPM for all other sizes	<ul style="list-style-type: none"> • Completion of finish • Removal of scratches • Use flex and pressure 	<ul style="list-style-type: none"> • Available Regular center sizes: 1/2", 5/8", & 3/4" • Available Mini center sizes: 3/8", 1/2", & 5/8"
Fine FlexiDisc® 	10,000 to 15,000 RPM for all sizes	<ul style="list-style-type: none"> • Beginning final polish • Use flex and pressure 	<ul style="list-style-type: none"> • Available Regular center sizes: 1/2", 5/8", & 3/4" • Available Mini center sizes: 3/8", 1/2", & 5/8"
Superfine FlexiDisc® 	10,000 to 15,000 RPM for all sizes	<ul style="list-style-type: none"> • Enhancing final polish • Use flex and pressure 	<ul style="list-style-type: none"> • Available Regular center sizes: 1/2", 5/8", & 3/4" • Available Mini center sizes: 3/8", 1/2", & 5/8"
NEW Medium FlexiCups™+ & FlexiPoints™+ 	3,000 to 8,000 RPM with Light Pressure for all cups and points • Do NOT exceed 8,000 RPM	<ul style="list-style-type: none"> • Use with water • Only use from top, not side 	<ul style="list-style-type: none"> • Achieve labial characterization
NEW Superfine FlexiCups™+ & FlexiPoints™+ 	3,000 to 8,000 RPM with Light Pressure for all cups and points • Do NOT exceed 8,000 RPM	<ul style="list-style-type: none"> • Use with water • Only use from top, not side 	<ul style="list-style-type: none"> • Polish gingival margin and labial surfaces
FlexiBuff® 	15,000 RPM for all sizes	<ul style="list-style-type: none"> • Add Enamelize® to tooth first and use pressure 	<ul style="list-style-type: none"> • Available Regular center sizes: 5/8" • Available Mini center sizes: 3/8", 1/2", & 5/8"