EN

**COSMEDENT** www.cosmedent.com | 800.621.6729  
**GINGA-Guard®** | LIGHT-CURE RESIN DAM

## INSTRUCTIONS FOR USE

## DESCRIPTION

**GINGA-Guard®** is a light-cure resin used to protect the gingival tissue from irritating materials such as hydrogen peroxide. It can be used as an alternative to a rubber dam for protecting the gingival surface.

**GINGA-Guard** cures with an LED or Halogen light at wavelength range of 400 - 500 nm.

## COMPOSITION

**GINGA-Guard** consists of Acrylic Polymer (73%). It is filled with inorganic filler and silicon dioxide (26%). Initiators and pigments are additional ingredients (<1%)

## INDICATIONS

- Protecting the gingival tissue during power bleaching.
- Protecting adjacent tooth structure or the surrounding restorations during etching with microetchants or phosphoric acid.

## DISPENSING

Place the disposable tip over the syringe opening and screw it on firmly before dispensing the material. Throw away tip after use and immediately replace the cap on the syringe to avoid prepolymerization.

## DIRECTIONS FOR USE

- 1 Apply **GINGA-Guard** to areas you wish to protect, i.e. gingival areas, adjacent tooth structure or surrounding restorations. Make sure the gingiva is dry before applying **GINGA-Guard**.
- 2 After application, polymerize for 5 seconds on each surface.
- 3 To remove **GINGA-Guard**, gently lift the material with a spoon excavator; it will come away in one piece.

## HYGIENE/DISINFECTION

- Disinfection of the syringe can be carried out using commercially available disinfecting solution.
- All instruments used should be thoroughly washed and autoclaved.

## IMPORTANT PRODUCT PRECAUTIONS

- **GINGA-Guard** tips are disposable and should only be used for one patient. Cover **GINGA-Guard** syringe with a cap immediately after use.
- Disinfect **GINGA-Guard** syringe between each patient.

## STORAGE

- Store at 39-75° F (4 - 24°C).
- Shelf life: see information on syringes and/or packages.
- Keep away from sunlight.
- Keep material out of reach of children.
- Do not use after the expiration date.

**Rx** IMPORTANT  
PRODUCT  
INFORMATION

Federal (U.S.) law restricts this device to sale by or on the order of a dentist.

**Note:** This product was developed for use in the dental sector and should be used in accordance with the instructions for use. The manufacturer does not accept liability for damage caused by its use for any other purpose. The user is also obliged to accept responsibility for checking before use that the material is suitable and can be used for the intended purpose, particularly if this is not listed in the instructions for use.

## DISPOSAL INSTRUCTIONS

Dispose of contents/ container in accordance with local regulation.

## CONTRAINDICATIONS

**GINGA-Guard** is contraindicated if patient is known to be allergic to any of the ingredients of **GINGA-Guard**. It is also contraindicated if the gingiva is inflamed.

## SIDE EFFECTS

Localized overheating may occur from overexposure to the light curing lamp.

## INTERACTIONS

None

## NOTICE

During time-consuming procedures, avoid overexposure from the operatory light to protect against premature polymerization of the material. The light curing unit must be maintained and calibrated according to the manufacturer's specifications. Serious incidents must be reported to DeltaMed GmbH, the responsible importer and the responsible authority.

## WARNING

**GINGA-Guard** contains Urethane acrylate and may cause irritation or sensitization in some people. In case of contact with the eyes, rinse immediately with plenty of water and seek medical attention.

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**Cosmedent® Inc.**