

GLOSS™ Light-Cure Temporary Glaze

DIRECTIONS FOR USE

PRODUCT DIRECTIONS:

Gloss™ is a light-cure one component varnish that does not leave a smear layer when cured. **Gloss** is used for coating temporary crown and bridge materials based on acrylic systems.

Properties

The special properties of **Gloss** include a high initial gloss, a full dry surface, a non-existent smear layer after light polymerization, and a more abrasion resistant surface.

Gloss is resistant to the adhesion of food and plaque.

APPLICATION

1. Remove the smear layer from your temporary crown and bridge material with a suitable solvent e.g. alcohol.
2. Apply a thin layer of **Gloss** with a brush onto the dull and dry surface, making sure there are no air bubbles.
3. Light-cure **Gloss** immediately after application from each side: Dental-Halogen Light Cure Unit (min. 500 mW/cm²): 20 sec
LED Curing Light (min. 1000 mW/cm²): 20 sec.
4. If required, after light curing clean the surface with alcohol and water. The surface will shine to a very high gloss.

Alternatively **Gloss** coating can be cured in a light oven for at least 90 seconds.

Note: Close bottle immediately after use.

Contraindications

Gloss is free of methylmethacrylate.

Do not use with patients with known allergy to methacrylate. In case of an allergic reaction please contact a physician.

Basic Safety Instructions

- **Gloss** is flammable. Avoid sources of ignition and do not store in direct sunlight.
- In case of contact with oral tissue or skin, flush immediately with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Keep away from children.
- **Caution: Federal law restricts the sale of this device by or on the order of a dentist.**

Storage

- Store at 40° – 75° F (4° – 24° C).
- Do not use after expiration date.
- Close bottle immediately after use to avoid exposure to light.

GLOSS™
Light-Cure Temporary Glaze
WWW.COSMEDENT.COM 800.621.6729