Insure® / Insure® Lite / Insure® Color Modifier | UNIVERSAL CEMENTATION SYSTEM

INSTRUCTIONS FOR USE

DESCRIPTION

Insure® is a dual cured, microhybrid filled, radiopaque, dental resin cement according to EN ISO 4049 (Type2/Class2/Group1 and Type2/Class3). Insure resin cement is available in two viscosities. Insure resin cement in the standard kit is formulated in a medium, thixotropic viscosity. Insure Lite resin cement is produced in a lower, more flowable viscosity. It is recommended that Insure be used for ceramic laminate veneers and Insure Lite be used for ceramic or composite crowns, inlays, onlays, and bridges. Insure cures with light at wavelength range of 400 – 500 nm.

Insure Simulcure Catalyst can be mixed with the colored resin cement in equal parts to produce a complete chemical cure in areas that are less accessible to light-curing. If using the dual-cure catalyst with Insure / Insure Lite the working time is 3 – 4 minutes and setting time is 5 – 6 minutes. Insure Lite Clear is available in an Automix version.

SHADES

Insure and Insure Lite are available in the following shades:

Clear, Yellow Red Light, Yellow Red Universal, Grey Light, White Opaque, Yellow Opaque, and Pink Opaque.

RELIABLE & PREDICTABLE

Insure resin cements have over 15 years of reliable performance and top ratings by the major dental materials advisories. This versatile resin cement system allows you to cement in veneers, inlays, onlays, crowns, metal and posts and cores with complete predictability. Insure features easy setting and cleanup, a choice of two viscosities, and a wide color range for predictable shade adjustments. Insure can be made dual-cure for inlays, onlays and crowns with the addition of Insure Simulcure Catalyst / Insure Lite Simulcure Catalyst.

COMPOSITION

Insure contains Bis-GMA, UDMA and 1,4-Butanediol dimethacrylate (30%). The fillers contain barium glass and silicon dioxide (69% by weight and 46% by volume). Initiators and stabilizers are additional ingredients (< 1%). The particle size of the fillers range from 0.04 – 3 micron.

Insure Lite and Insure Color Modifier contain Bis-GMA, UDMA and 1,4-Butanediol dimethacrylate (35%). The fillers contain barium glass and silicon dioxide (64% by weight and 41% by volume). Initiators and stabilizers are additional ingredients (< 1%). The particle size of the fillers range from 0.04 – 3 micron.

Insure Simulcure Catalyst / Insure Lite Simulcure Catalyst contain Bis-GMA, UDMA and 1,4-Butanediol dimethacrylate (30% / 35%). The fillers contain barium glass and silicon dioxide (69% / 64% by weight and 46% / 42% by volume). Both contain benzoyl peroxide. The particle size of the fillers range from 0.04 – 3 micron.

INDICATIONS

(Areas of Application)
– Cementation of ceramic veneers
– Cementation of ceramic or composite crowns, inlays, onlays, and bridges

DISPENSING INSURE

Place the disposable tip over the syringe opening and screw it on firmly before dispensing the Insure resin cement. This will provide controlled placement into the veneer shell without excess resin cement extruding from the syringe causing material waste. Throw away tip after use and immediately replace the cap on the Insure syringe to avoid pre-polymerization.

EVALUATING FOR MARGINAL FIT

Prior to cementation of ceramic laminate veneers, it is advisable to evaluate each ceramic laminate veneer for precise marginal fit and interproximal adaptation to adjacent ceramic laminate veneers or natural tooth structure. Any aberrations should be identified, and appropriate measures set up for their elimination.

SELECTING SHADE

We recommend that you use Prevue® Try-In, Cosmedent's water soluble matching try-in gels, to select the correct shade of your Insure resin cement. Prevue Try-In Gels are especially effective for feldspathic porcelain veneers or any minimal or no-prep veneers. These non-light activated pastes accurately match the polymerized color of our seven resin cements and wash away quickly with water, making shade selection predictable and fast.

PLACEMENT DIRECTIONS

Insure is recommended for ceramic veneers and Insure Lite (thin viscosity) is recommended for ceramic or composite crowns, inlays, onlays, bridges and veneers because of the optimum film thickness.

Once the ceramic laminate veneers have been tried in and approved for marginal fit, proximal adaptation, and resin cement color, the following cementation procedure is undertaken:

- 1 Isolate the teeth to be veneered with a rubber dam or cotton rolls and an effective evacuation system that will sustain a dry uncontaminated field throughout the duration of the cementation procedure.
- 2 Place retraction cord, if needed, along the cervical margins of all teeth to be veneered to permit clear visualization of the cervical finish lines and counteract the potential for soft tissue contamination of any kind.
- 3 Floss each contact to remove any residual provisional veneer resin or resin cement. Relax but do not open the contacts with controlled passage of narrow, Superfine, metal backed FlexiStrips™ (Cosmedent). Remove any residual provisional veneer cement on the facial surface with a fine diamond bur.
- 4 Clean all tooth surfaces to be veneered by initially pumicing all accessible tooth surfaces with pumice. Rinse away all pumice residue and microetch any existing resin restorations to which the veneers will be cemented. Rinse all tooth surfaces to be cemented thoroughly with water and air dry.
- 5 If the laboratory has not acid etched and silanated the restorations or if this surface is altered by the dentist, then these steps must be performed prior to resin bonding. If the laboratory has etched the porcelain, skip to step 5c.
 - a) Etch the intaglio of the veneer for 20 seconds with the 5% Porcelain etchant (e.max, Ivoclar Vivadent) or 3 minutes with 9.5% Hydrofluoric acid if Feldspathic porcelain.
 - b) Rinse thoroughly with water and air dry with oil free air to visualize a frosted, etched appearance.
 - c) Apply a silanating agent such as Insure Silanator (Cosmedent) for 30 seconds to all etched inner surfaces of the veneers and gently air dry. Place the ceramic laminate veneers in sequential order for cementation in pairs, starting with the central incisors, or in quartets, starting with the incisors.
- Place a thin metal strip interproximally on the distal proximal surface of the last prepared tooth to be cemented in this grouping (either a pair of veneers or 4 consecutive veneers starting at the midline), to prevent etching an adjacent tooth not in this cementation grouping. Etch each tooth surface to be veneered with 35%- 37% phosphoric acid for approximately 20 seconds. Rinse off all acid with a gentle water spray for 5 seconds and carefully dry the teeth with a clean air spray to a damp, but not desiccated state.

or totally set up before subsequent veneers are ready to be placed. Working time for Insure / Insure Lite is 3 – 4 minutes and setting time is 5 – 6 minutes.

HYGIENE/DISINFECTION

– Disinfection of the syringe can be carried out using a commercially available disinfecting solution.
– All instruments used should be thoroughly washed and autoclaved.

STORAGE

– Once opened, do not refrigerate
– Keep Insure at room temperature (approximately 39 – 75°F; 4 – 24°C)
– Close syringe immediately after use. Exposure to light causes premature polymerization
– Do not use after expiration date
– Keep away from sunlight
– Keep material out of reach of children.

IMPORTANT PRODUCT INFORMATION

Federal (U.S.) law restricts this device to sale by or on the order of a dentist.

Note: This product was developed for use in the dental sector and should be used in accordance with the instructions for use. The manufacturer does not accept liability for damage caused by its use for any other purpose. The user is also obliged to accept responsibility for checking before use that the material is suitable and can be used for the intended purpose, particularly if this is not listed in the instructions for use.

DISPOSAL INSTRUCTIONS

Dispose of contents/ container in accordance with local regulation.

CONTRAINDICATIONS

Insure is contraindicated if a patient is known to be allergic to any of the ingredients of Insure.

SIDE EFFECTS

May cause irritation in sensitive patients.

INTERACTIONS

Do not use in combination with eugenol containing materials as these may impede the polymerization process.

NOTICE

During time-consuming procedures, avoid overexposure from the operatory light to protect against premature polymerization of the material. The light curing unit must be maintained and calibrated according to the manufacturer's specifications. Serious incidents must be reported to DeltaMed GmbH, the responsible importer and the responsible authority.

WARNING

Insure contains 1,4-Butanediol dimethacrylate, Urethane dimethacrylate and Di-benzoyl peroxide (only in Simulcure Catalyst). May cause an allergic skin reaction. Causes skin irritation. Causes serious eye irritation. Wear protective gloves. In case of contact with the eyes, rinse immediately with plenty of water and seek medical attention.

PACKAGING

Insure 3 g syringes (includes 3 Nozzle Tips White)
Insure Light 3g syringes (includes 3-17 gauge tips)
Insure Color Modifier (includes 3-19 gauge tips)

Canada Site Licence Holder:

CLINICAL RESEARCH DENTAL
167 Central Ave, Suite 200
London, Ontario, Canada N6A 1M6
Establishment License No. 666

Australia Sponsor:

ERSKINE ORAL CARE
5-9 Kylie St,
Macksville NSW 2447
Australia

Australia Sponsor:

Amalgadent Dental Supplies Pty Ltd
334 Glen Huntly Road
Elsterwick, Victoria 3185 Australia

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