EN

COSMEDENT™ www.cosmedent.com | 800.621.6729

Oxygone™

INSTRUCTIONS FOR USE

DESCRIPTION

Oxygone™ is a water-soluble gel used to prevent the air inhibited layer from forming on cured composites and luting cements.

COMPOSITION

Oxygone consists of propylene glycol, glycerin, water and acrylic polymer. It is pH neutral.

INDICATIONS FOR USE

(Areas of application)

- To cover composite filling materials and luting cements during curing to prevent the formation of an air inhibited layer.

DISPENSING OXYGONE

Place the disposable tip over the syringe opening and screw it on firmly before dispensing **Oxygone**. Throw away tip after use and replace the cap on the **Oxygone** syringe.

APPLICATION

DIRECT COMPOSITE

After modelling the desired anatomical contour of the restoration, coat the entire surface to be cured with **Oxygone**. Light cure the composite for the final cure according to the manufacturer's instructions. After the final cure, rinse **Oxygone** from the surface with a water spray.

PORCELAIN VENEER CEMENTATION

After seating the veneers, tack cure the cement. Coat the entire marginal surface of the restoration-interface with **Oxygone** and perform the final cure. Rinse **Oxygone** from the surface with a water spray.

HYGIENE/DISINFECTION

- Disinfection of the syringe can be carried out using a commercially available disinfecting solution.
- All Instruments used should be thoroughly washed and autoclaved.

STORAGE

- Store at 39 – 75°F (4 – 24°C).
- Shelf life: see information on syringes and /or packages.
- Keep away from sunlight
- Keep material out of reach of children.
- Do not use after the expiration date.

Rx IMPORTANT
PRODUCT
INFORMATION

Federal (U.S.) law restricts this device to sale by or on the order of a dentist.

Note: This product was developed for use in the dental sector and should be used in accordance with the instructions for use. The manufacturer does not accept liability for damage caused by its use for any other purpose. The user is also obliged to accept responsibility for checking before use that the material is suitable and can be used for the intended purpose, particularly if this is not listed in the instructions for use.

DISPOSAL INSTRUCTIONS

Dispose of contents/ container in accordance with local regulation.

CONTRAINDICATIONS

Placement of **Oxygone Gel** is contraindicated if a patient is known to be allergic to any of the ingredients.

SIDE EFFECTS

May cause irritation from chronic exposure.

INTERACTIONS

None

NOTICE

Serious incidents must be reported to DeltaMed GmbH, the responsible importer and the responsible authority.

WARNING

Oxygone Gel contains glycerin and propylene glycol and may cause irritation in sensitive persons. In case of contact with the eyes, rinse immediately with plenty of water and seek medical attention if irritation persists.

PACKAGING

Oxygone Gel is available in 3 g syringes (includes 3-25 gauge tips).

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Schwerwiegende Vorfälle sind der DeltaMed GmbH, dem zuständigen Importeur und der zuständigen Behörde zu melden.
Cosmedent® Inc.