EN

COSMEDENT® www.cosmedent.com | 800.621.6729

RSVP™

INSTRUCTIONS FOR USE

DESCRIPTION

RSVP™ (Rapid Simplified Veneer Provisionals) is an innovative restoration system used to create veneer provisionals. **The RSVP System** is composed of two specially formulated light-cure microhybrid composite resins: **RSVP Lite** and **RSVP Heavy**. **RSVP Lite** is a flowable composite and **RSVP Heavy** is a sculptable stiff composite. Both cure with light at wavelength range of 400 – 500 nm.

COMPOSITION

RSVP Lite contains UDMA and 1,4-Butanediol dimethacrylate (46%). The fillers contain strontium glass, silicon dioxide and acrylic polymer (53% by weight and 37% by volume). Initiators, stabilizers and pigments are additional ingredients (< 1%).

RSVP Heavy contains UDMA, Bis-GMA and 1,4-Butanediol dimethacrylate (31%). The fillers contain strontium glass, silicon dioxide and acrylic polymer (68% by weight and 51% by volume). Initiators, stabilizers and pigments are additional ingredients (< 1%).

INDICATIONS FOR USE

(Areas of application)
– Temporary or provisional veneers

DISPENSING COMPOSITE MATERIAL**Heavy Viscosity:**

– Extrude composite onto a mixing pad by turning the handle slowly in a clockwise manner. Slice an increment of proper size with a sterile instrument onto a mixing pad.
– Immediately replace syringe cap to prevent polymerization. If the extruded material is not used right away, cover the composite with a light protecting box such as our **ResinKeeper™**.

Lite Viscosity

– Place the disposable tip over the syringe opening and screw it on firmly before dispensing the **RSVP Lite**.
– Throw away tip after use and immediately replace the cap on the **RSVP** syringe to avoid prepolymerization.

APPLICATION

For a detailed description and technique video please visit Cosmedent.com.

- Take an **RSVP** polysiloxane impression of a pre-treatment model using a transparent, non-perforated full arch tray to register all teeth to be provisionalized plus two posterior to the last tooth on each side.
- Remove the impression tray from the model; the impression remains on the model. Trim the impression using a #15 scalpel 1 mm from the papilla tips on the facial and 3 mm from the incisal edges on the palatal (occlusal).
- Remove the impression from the model and seat it in the patient's mouth.
- Isolate the teeth to be provisionalized with cotton rolls and dry. The impression is filled with **RSVP Lite**. It is advisable to place a light coat of an unfilled light-cure resin on the dry tooth surfaces and light cure each tooth 5 – 10 seconds.
- Seat the impression and have an assistant exert light pressure from the occlusal.
- Remove excess **RSVP** with a fine brush. **RSVP** should reach the rims of the impression. Light cure each tooth from the facial for 5 seconds followed by curing for 5 seconds from the incisal.
- Remove the impression from the provisionals. Extend the cervical border with additional **RSVP Lite** to no closer than 1 mm to the cervical margin. Light cure the extension for 10 seconds per tooth.
- Place **RSVP Heavy** in the form of a flat patty directly on the middle of the cervical one-third of the provisional. Feather cervically and incisally with a thin composite instrument to approximately ½ mm overcontoured on the facial. Light cure 10 seconds.
- Finish and polish in the usual manner.

HYGIENE / DISINFECTION

– Disinfection of the syringe can be carried out using a commercially available disinfecting solution.
– All instruments should be thoroughly washed and autoclaved.

STORAGE

– Store at 40 – 75°F (4 – 24°C).
– Shelf life: see information on syringes and / or packages.
– Keep away from sunlight.
– Keep out of reach of children.
– Do not use after the expiration date.

IMPORTANT PRODUCT INFORMATION

Federal (U.S.) law restricts this device to sale by or on the order of a dentist.

Note: This product was developed for use in the dental sector and should be used in accordance with the instructions for use. The manufacturer does not accept liability for damage caused by its use for any other purpose. The user is also obliged to accept responsibility for checking before use that the material is suitable and can be used for the intended purpose, particularly if this is not listed in the instructions for use.

DISPOSAL INSTRUCTIONS

Dispose of contents/ container in accordance with local regulation.

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Cosmedent® Inc.
DeltaMed GmbH