EN

**COSMEDENT**  
www.cosmedent.com | 800.621.6729

## COSMECORE®

DUAL-CURE AUTOMIX CORE BUILD-UP MATERIAL

### INSTRUCTIONS FOR USE

#### GENERAL CHARACTERISTICS

COSMECORE® is a dual-cure, automix, hybrid filled, radiopaque, composite core buildup material according to EN ISO 4049 (Type 2/Class 3). The automix syringe is a dual cartridge system. The base and catalyst are automixed by forcing the materials through a disposable mixing tip. COSMECORE is light cured at wavelength range of 400 – 500 nm and self-cures in the absence of light.

**Please note:** Although COSMECORE is dual-cure, it must be light-cured to ensure quality.

#### COMPOSITION

The monomer matrix of COSMECORE contains UDMA, Bis-EMA, Bis-GMA, Tricyclodecanedimethanol diacrylate and 1,4-Butanediol dimethacrylate (30%). The fillers contain barium glass and silicon dioxide (69% by weight and 48% by volume). Initiators, stabilizers and pigments are additional ingredients (< 1%). The particle size of the fillers range from 0.04 – 5 micron. The catalyst contains Dibenzoyl peroxide.

#### TIMING

##### SELF-CURING MODE

**Working time:** 90 Seconds

**Curing time:** 4 – 5 Minutes

#### INDICATIONS FOR USE

– Core buildups

#### DISPENSING COSMECORE

It is **very important** to clear or bleed the cartridge or syringe **prior to applying the mixing tip**. Gently press the dispenser lever until both components (base and catalyst) begin to flow out evenly onto a mixing pad.

To apply the mixing tip, push the tip completely down until the notch of the mixing tip is aligned with the notch on the housing of the cartridge/syringe. Secure the mixing tip in place by gripping the colored base, (not the clear mixing tip), and by turning it 1/4 turn clockwise.

**For each application, place a new mixing tip on the cartridge/syringe.** For each new tip, bleed a small pea size amount of material prior to use.

#### PLACEMENT DIRECTIONS

① Etch the tooth. Rinse, but leave the tooth moist.

② Rinse, but leave the tooth moist.

③ For maximum results and a stronger bond, apply an advanced generation bonding agent such as dual-cure **Complete™** or a dual-cure 4th generation bonding agent, such as ALL-Bond 3® (Bisco, Inc.). **Because of the acidic nature of most 5th generation bonding agents, using them may reduce bond strength and may cause your buildups to debond.**

④ Dispense the material from the automix cartridge into the preparation and **light cure for 5 seconds. WAIT 5 seconds**, then **continue to light cure for an additional 40 seconds. This procedure will reduce heat buildup.**

⑤ The buildup is ready to prep after light curing.\*  
\*The dual-cure capabilities ensure that the material will continue to self-cure after light curing.

#### HYGIENE / DISINFECTION

To prevent polymerization or clogging, the dual-syringe mixing tip should be left in place until the next application. If used intraorally, wipe mixing tip with disinfectant prior to storage. The Automix-Dispenser is reusable after disinfection. Disinfection can be carried out using a commercially available disinfecting solution.

#### STORAGE

– Once opened, do not refrigerate.  
– Store at room temperature, approximately 39 – 75°F (4 – 24°C).  
– Shelf life: see information on syringes and/or packaging.  
– Keep away from sunlight.  
– Keep out of reach of children.  
– Do not use after the expiration date.

#### **Rx** IMPORTANT PRODUCT INFORMATION

Federal (U.S.) law restricts this device to sale by or on the order of a dentist.

**Note:** This product was developed for use in the dental sector and should be used in accordance with the instructions for use. The manufacturer does not accept liability for damage caused by its use

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**Cosmedent® Inc.**
**DeltaMed GmbH**