EN

COSMEDENT™ www.cosmedent.com | 800.621.6729

Renamel® De-Mark™ FLOWABLE-HYBRID RADIOPAQUE LINING MATERIAL

## INSTRUCTIONS FOR USE

## DESCRIPTION

**Renamel® De-Mark™** is a hyper-radiopaque flowable cavity liner that clearly demarcates your composite restoration on an x-ray. It flows easily and adapts perfectly to the floor of the prep and the gingival seat areas. The use of **Renamel De-Mark** will help the practitioner easily diagnose recurrent decay.

**Renamel De-Mark** is a light-cured, microhybrid filled, radiopaque (220% Al) dental filling material, according to EN ISO 4049 (Type 2/Class 2/ Group 1). **Renamel De-Mark** cures with an LED or Halogen light at wavelength range of 400-500 nm.

## COMPOSITION

**Renamel De-Mark** contains Bis-GMA, UDMA and 1,4-Butanediol dimethacrylate (33%). The fillers contain barium glass, silicon dioxide and barium sulfate (66% by weight and 42% by volume). Initiators and stabilizers are additional ingredients (< 1%). The particle size of the fillers range from 0.04 – 7 micron.

## INDICATION FOR USE

- Cavity liner
- Radiopaque lining material to help evaluate marginal integrity, rule out overextended restorations, and to help detect recurrent caries.

## DISPENSING

- Place the disposable tip over the syringe opening and place it on firmly before dispensing the material.
- Throw away tip after use and immediately replace the cap on the syringe to avoid prepolymerization.

## DIRECTIONS FOR USE

- Following polymerization of bonding agent, place a thin layer of **Renamel De-Mark** at that gingival seat area and polymerize 40-60 seconds. This reduces polymerization shrinkage and postoperative sensitivity.
- Place restorative material of choice, cure, finish and polish. Note: For deep cavities with only a minimal layer of dentin over the pulp or those which have already had a small area exposed, the pulp should be indirectly or directly capped locally using calcium hydroxide compound. The calcium hydroxide compound should be covered with a glass ionomer cement before the **Renamel De-Mark** is inserted into the cavity.

## HYGIENE/DISINFECTION

- Disinfection of the syringe can be carried out using commercially available disinfecting solution.
- All instruments used should be thoroughly washed and autoclaved.

## STORAGE

- Once opened, do not refrigerate.
- Store at 39 – 75°F (4 – 24°C).
- Shelf life: see information on syringes and/or packaging.
- Keep away from sunlight.
- Keep out of reach of children.
- Do not use after the expiration date.

## Rx IMPORTANT PRODUCT INFORMATION

Federal (U.S.) law restricts this device to sale by or on the order of a dentist.

**Note:** This product was developed for use in the dental sector and should be used in accordance with the instructions for use. The manufacturer does not accept liability for damage caused by its use for any other purpose. The user is also obliged to accept responsibility for checking before use that the material is suitable and can be used for the intended purpose, particularly if this is not listed in the instructions for use.

## CONTRAINDICATIONS

**Renamel De-Mark** is contraindicated if a patient is known to be allergic to any of the ingredients in **Renamel De-Mark**.

## SIDE EFFECTS

May cause irritation in sensitive patients.

## INTERACTIONS

Do not use with eugenol containing materials, as these may impede the polymerization process.

## NOTICE

During time-consuming procedures, avoid overexposure with the operatory light to protect against premature polymerization of the material. The light curing unit must be maintained and calibrated according to the manufacturer's specifications. Serious incidents must be reported to DeltaMed GmbH, the responsible importer and the responsible authority.

## DISPOSAL INSTRUCTIONS

Dispose of contents/ container in accordance with local regulation.

## WARNING

**Renamel De-Mark** contains 1,4-Butanediol dimethacrylate and Urethane dimethacrylate. Causes skin irritation. Causes serious eye irritation. May cause an allergic skin reaction. Wear protective gloves. In case of contact with the eyes, rinse immediately with plenty of water and seek medical attention.

## PACKAGING

- **Renamel De-Mark** 3 gram Syringe – #7011DM (includes 3-18-gauge and 3-19-gauge tips)
- **Renamel De-Mark** Refill Tips (18 gauge) Package of 20 – #203102
- **Renamel De-Mark** Refill Tips (19 gauge) Package of 20 – #203103

Canada Site Licence Holder:  
CLINICAL RESEARCH DENTAL  
167 Central Ave, Suite 200  
London, Ontario, Canada N6A 1M6  
Establishment License No. 666

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**Cosmedent® Inc.**